

The Aeria 5+ is our premium mattress overlay system designed to meet the acute care needs of home or long term care environments.

Indicated for the prevention and treatment of skin breakdown and

pressure ulcers in patients at medium risk, the overlay is constructed from 18 transverse cells – 5 inch 120 mm high – that cyclically inflate and deflate in an alternating pattern.

Cyclic alternation of pressure prevents arterial and venous capillary occlusion in surface tissue – maintaining and simulating the flow of blood and lymphatic fluids through these tissues to provide essential oxygen and remove metabolic waste.



control unit

- Digital control unit pumps over six liters of air per minute; easy to operate digital touch panel with bright LED indicators.
- Adjustable pressure settings to match individual user weight and comfort needs.
- Quick connector air hose attachments.
- Visible and audible alarm warns of power failure or when cell pressure drops below the accepted minimum level.
- Reinforced air hoses prevent unwanted kinks or interference to air flow.
- Inbuilt rear hanging hooks for easy mounting and control unit portability.

custom designed carry bag

Stylish carry bag for compact transport or storage. Separate compartments store all product components and collateral, including simple usage instructions and our easy to follow Quick Set Up Guide.

overlay

- Vapor permeable, waterproof top cover keeps moisture away from the skin; multi stretch fabric minimizes friction and reduces shear.
- Seams are sewn and tape sealed for added bacterial protection.
 - 360 degree, fully enclosable zipper with quick release functionality makes it easy to remove the top cover from a standing position.
- 18 x 5 inch 120 mm high density cells made from premium quality TPU coated nylon for superior performance, wear and tear. Each cell connects independently to mattress base for easy removal for cleaning or replacement.
- Rapid release CPR tag for emergency deflation.
 - Static head cells for stable head support and optimum user comfort.
- Static transport function for stable patient handling and transfer.
- Strong elastic straps simply and securely lock the overlay to any standard single bed mattress.

technical specification 220V

Your local
Aeda Healthcare
authorized dealer:

general	product type	dynamic mattress overlay
	model no	Aeria 5+ (MT-A1-151)
	alternation sequence	1 in 2
	alternation cycle	10 ~ 12 minutes
	maximum user weight	350 lb 160 kg
	CPR function	mattress tag
	low air loss feature	✗
	transport function	✓
	static mode	✓
	risk category	medium
	carry bag	✓
	certifications	CE, UL
	warranty	2 years
	contraindications	instable spinal cord injury
cervical traction		
intended care settings	long term or extended care	
	home care	
control unit	style	digital
	length	11.5" 290 mm
	width	5.5" 140 mm
	height	10.5" 270 mm
	weight	6.5 lb 2.9 kg
	air output	6 liters per minute
	pressure range	21 ~ 59 mmHg
	alarm mode	visible / audible
	quick release connectors	✓
	rated voltage	AC110V 60Hz 220V 50Hz
	rated input power	20 va
	protection class	Class 2
	power cord length	16 ft 5 m
	mattress	length
width		35" 880 mm
height		5" 120 mm
weight		11.5 lb 5.3 kg
top cover material		PU coated multi-stretch polyester
skirt seams		sewn & tape sealed
cover attachment		quick release zipper
base material		PU coated nylon
non skid pad		✗
lifting handle		✓
anti-kink hosing		✓
cells		no. of cells
	cell height	5" 120 mm
	cell material	TPU coated nylon
	cell in cell	✗
	ventilated air cells	✗
	static head cells	✓

Sep09